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**BRESSER Science IVM-401
Transmission-type Inverted Microscope**

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Caution!

To work with this microscope, sharp and pointed aids are being used. Please take care that this microscope and its accessories are stored at a place out of reach of children. Let children only work with this microscope under an adult's supervision! Keep packing material (plastic bags etc.) away from children!

**Do not dispose of electric equipment together with household waste material!**


In observance of European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE) and its implementation in accordance with national law, electric equipment that have reached the end of its life must be collected separately and returned to an environmentally compatible recycling facility.

Discharged batteries and damaged re-chargable batteries must be disposed of at special battery collection points. Information is available from your local disposal agent or local authority regarding the disposal of devices or batteries manufactured after the 01.06.2006.

Note!

If you have any complaints or queries please first contact your national service centre by telephone. The address is included in these instructions.

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Telefon: _____

Kaufdatum: _____

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The model "BRESSER Science IVM-401" Inverted Microscope has the feature that observations in a culture bottle or a culture utensil can be done. One can observe transparent and protozoan living specimens. It is an advanced instrument, suitable for the research of living specimen structure, liquids,

deposits and so on. It also can be used in research institutes, universities, medical treatment, agriculture and animal husbandry etc.

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- I. Components
- II. Technical Data
- III. Observing Operation
- IV. The Operations of Each Unit
- V. Maintenance, EEC Conformity Explanation, Warranty

- 14 Condenser adjustment screw
- 15 Condenser stage tray
- 16 Aperture diaphragm adjustment swivel
- 17 Lamp holder

I. COMPONENTS

Fig. 1

- 1 Lamp box
- 2 Filter seat
- 3 Annular diaphragm slab
- 4 Annular diaphragm adjustment screw
- 5 Eyepiece
- 6 Trinocular head
- 7 Knob for tensional adjustment of focusing
- 8 Coarse focusing knob
- 9 Power supply plug
- 10 Mains in
- 11 Fine focusing knob
- 12 Objective
- 13 Stage

Abb. 2

- 18 Lamp holder adjustment pole
- 19 Collector adjustment lever
- 20 Field diaphragm adjustment lever
- 21 Condenser holding screw
- 22 Reticle
- 23 Condenser up/down knob
- 24 Stage moving system
- 25 Longitudinal adjustment knob
- 26 Transversal adjustment knob
- 27 Power switch
- 28 Brightness control knob
- 29 Holding screw for stage moving system
- 30 Limit knob (up stop)
- 31 Tube holding screw
- 32 Rod for observation/photo changing
- 33 Photo tube

II. TECHNICAL DATA

Eyepieces

Eye-piece type	Magnification	Field of view \varnothing / mm	Focal length f / mm	Remarks
Wide field eyepiece (WF)	10x	20	25	
WF-eyepiece with crosshair	10x			optional

Objectives

Objective type	Magnification	Numerical aperture NA	Working distance WD / mm	Coverglass thickness d / mm	Remarks
Long working distance planachromatic objective	10x	0.25	8.1	-	
	25x	0.40	4.8	1.2	
	40x	0.60	3.3	1.2	
Long working distance plan phase contrast objective	10x	0.25	8.1	1.2	optional erhältlich
	25x	0.40	4.8	1.2	
	40x	0.60	3.3	1.2	optional erhältlich

Total magnification

Objectives:	10x	25x	40x
Total magnification			
Eye-piece:			
10x	100x	250x	400x

- Condenser (with phase contrast): Working distance 30 mm.
- Stage: Moveable range 79 mm x 112 mm.
- Coaxial coarse/fine focusing system with limit knob and tensional adjustment ring: Minimum division of fine focusing scale 0.002 mm.

- Interpupillary distance adjustment range: 53 mm - 75 mm.
- Illumination system: 6 V, 30 W halogen bulb lamp with brightness adjustment.
- Power supply: 230 V (50 Hz).
- Anti-fungus.

III. OBSERVING OPERATION (Fig. 1 & 2)

1. Insert the condenser into the condenser stage tray (15) and fix it by the condenser holding screw (21). Plug the power supply plug (9) into the mains in (10) and the mains plug into the mains socket. Then push the power switch (27) onto "I". So, the power supply is switched on.
2. Put a specimen onto the center of stage. Turn the 10x objective to the working position, then adjust the coarse/fine focusing system.
3. Adjust the interpupillary distance and the diopter.
4. Adjust the position of the condenser, the brightness control knob (28) and the aperture diaphragm adjustment swivel (16) to get a satisfactory illumination state.
5. Little adjust the fine focusing knob (11), when you use objectives of different magnifications.

IV. THE OPERATIONS OF EACH UNIT

1. Adjustment of interpupillary distance (Fig. 3)


Put the specimen on the stage and bring it into the exact focus. Adjust the interpupillary distance of the binocular viewing by "folding" the binocular tube (2) on the optical head, until the right and the left field of view can be composed to one.

2. Adjustment of diopter (Fig. 3)


Put the specimen on the stage. Turn the 40x objective to the working position. Firstly, observe at the right eyepiece with the right eye; adjust the coarse and fine focusing knobs to image clearly. Secondly, observe at the left eyepiece with the left eye; adjust the diopter control (1) to image clearly.

3. Coarse/Fine focusing (Fig. 4)


The instrument uses a coaxial coarse/fine focusing mechanism. The knob for tensional adjustment of focusing (3) is used for adjusting the tension of the coarse focusing knob (2) to prevent the stage from naturally sliding down. The limit knob (4) prevents accidental contact between specimen and objective. After locking the limit knob (by turning up and fastening) at a chosen height of the stage, you cannot focus the object closer to the objective using the coarse focusing knob. So, the object is protected against damage. (1) is the fine focusing knob.

4. Stage (Fig. 5)

The seat of a culture utensil (4 or 5) can be put on the stage, used for different culture utensils and specimens. The longitudinal adjustment knob (1) and the transversal adjustment knob (2) of the mechanical stage are coaxial; adjust them very convenient. When you use a large culture utensil, you should loose the holding screw (3) for the stage moving system. Then take away the stage moving system from stage and put the large culture utensil on the stage directly.

5. Adjustment of lamp position (Fig. 6)

Turn the condenser up-down knob (Fig. 2, No. 23) to move the condenser to the reticle position (Fig. 2, No. 22), so that both marking lines are at the same height. Then put a piece of white paper on the filter seat (Fig. 1, No. 2). Adjust the collector adjustment lever (1), make the image of the filament clear on the paper. If the filament image is not on the center of the overpass hole, you should loosen the lamp holder fixing screw (2) and the lamp holder adjustment pole (Fig. 2, No. 18). Then move the lamp holder (3) in order to locate the image of the filament on the center of the overpass hole. Now, tighten the lamp holder fixing screw and the lamp holder adjustment pole to fix the lamp. So, you fix the position of the lamp.

6. Adjustment of the field diaphragm (Fig. 1 & 2)

Turn the 10x objective to the working position and observe with it. Turn the coarse/fine focusing knob (8 and 11) in order to make the image of the object clear. Turn the field diaphragm adjustment lever (20) to adjust the diminishing of the field diaphragm. Turn the condenser up-down knob (23) to gain the image of the field diaphragm, then adjust the condenser adjustment screw (14) to center the field diaphragm and the eyepiece diaphragm. Turn the field diaphragm adjustment lever (20) in order to fit the size of the field diaphragm to the eyepiece diaphragm.

7. Adjustment of condenser aperture diaphragm (Fig. 1)

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8. Power switch and brightness control (Fig. 2)


Push the power switch (27) to "I" position; now, the power is switched on. Then adjust the brightness control knob (28) to observe the image of the object comfortably.

NOTICE: Don't let the brightness control knob on the highest position for a long time, to avoid depressing the lifetime of the lamp.

9. Phase contrast attachment (Fig. 7 & 8)

1. Adjust the field diaphragm according to paragraph "6. Adjustment of the field diaphragm" and move the condenser to the corresponding position.
2. Turn the aperture diaphragm adjustment swivel (Fig. 1, No. 16) to open the aperture diaphragm maximally.
3. Turn the corresponding phase contrast objective (Fig. 1, No. 12) to the working position.
4. Take out one eyepiece, insert the centering telescope (1) into the eyepiece holder and adjust it. Make the image of the annular phase contrast dark halo (Fig. 8, No. 2) and of the annular diaphragm light halo (Fig. 8, No. 1) clear.


5. Adjust the annular diaphragm adjustment screw (Fig. 1, No. 4). Superposition the annular diaphragm light halo (Fig. 8, No. 1) and the annular phase contrast dark halo (Fig. 8, No. 2).
6. Take out the centering telescope and insert the eyepiece to observe with phase contrast.
7. If you change to a phase contrast objective of another magnification, you have to adjust according to above again.

10. Exchange of the lamp (Fig. 6)

1. Switch off the power switch (Fig. 2, No. 27) and pull out the mains plug.
2. Loosen the lamp holder fixing screw (2) and the lamp holder adjustment pole (Fig. 2, No. 18). Remove the lamp holder (3) from the lamp box.
3. Pull out the old lamp from the lamp holder (3). Be careful, as the lamp may be hot. Then insert the new lamp into the lamp holder, and clean the new lamp with a little alcohol. Let it dry.
4. The adjustment of the lamp position is done according to paragraph "5. Adjustment of lamp position".

11. Exchange of the fuse (Fig. 9)

Switch off the power switch (Fig. 2, No. 27) and pull out the mains plug from mains socket and the power supply plug (2) from mains in (3). Take out the fuse holder (1): Grip with a small screwdriver under the lug of the fuse holder (Fig. 9, arrow) and pull it carefully forwards. Change the fuse. Install the fuse holder and plug again.
 The specification of fuse: 250 V / 0,5 A, F.

12. Use of culture utensil (Fig. 1)

While you use the tall culture utensil for observation, push the condenser stage tray (15) from the light path.

V. MAINTENANCE

1. Sweep the lens

Sweep the lens by lens tissue or soft fabric immersed with a mixed liquid of alcohol/ether. Clean the 100x objective from oil whenever you finish operating.

2. Clean the painted parts

The dust on the painted parts can be removed by gauze. For the grease spots, the gauze immersed slightly with aviation gasoline is recommended. Do not use organic solvents such as alcohol, ether or other thinner etc. for cleaning the painted parts or plastic components.

3. Avoid disassembling the microscope

Because of being a precise instrument, do not disassemble the microscope casually. That may cause serious damage to its performance.

4. Being not used

Cover the microscope with the dust cover and place it there, where it is dry and mouldless. We suggest the storage of all objectives and eyepieces in a closed container with drying agent.

EEC CONFORMITY EXPLANATION

Bresser GmbH, resident in 46414 Rhede/Westf., Gutenbergstr. 2, Germany, explains the agreement with in the following specified EEC guidelines for this product:

EN 61326
 EN 2001+A2: 2005

Product description: Transmission-type Inverted Microscope
 Model: BRESSER Science IVM-401

Rhede, Juli 20, 2007

Bresser GmbH

Helmut Ebbert
 Managing director

WARRANTY

The period of warranty is 2 years, beginning on the day of purchase. Please keep the cash receipt as evidence of purchase. Devices which become defective during the warranty period can be returned to the dealer where the device was bought. The repaired device or a new one will then be returned to you. In the case of defects which occur after the end of the warranty period, the devices can also be returned. However, repairs which become necessary after the end of the warranty period will be subject to a service fee.

Important:

Make sure to return the device carefully packed in the original packaging in order to prevent transport damage. Please also enclose the cash receipt (or a copy). This warranty does not imply any restriction of your statutory rights.

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